Ezzi-Lift[®] PRODUCT SPECIFICATIONS

The Ezzi-lift microcurrent device for aesthetic use gently stimulates, tightens, and tones the face and body for a more youthful, natural appearance. The non-invasive, hand-held device can be used without creams or gels and is ideal for use at home or in the clinic.

BENEFITS

- Diminishes the appearance of fine lines and wrinkles
- Tightens skin around the neck and jawline
- Delivers gentle microcurrent massage to the face and body
- Enables treatment without the need for creams or gels

FEATURES

- Avazzia BEST[®] microcurrent platform
- Medical-grade stainless steel electrodes
- Multiple treatment modes and intensity options
- Accessory port for optional treatment tools

TECHNICAL SPECIFICATIONS

Size, weight	Size: 119 mm X 64 mm X 33 mm (4.69 in. x 2.52 in. x 1.30 in.) Weight (device only): Approx. 108 grams (3.8 oz) Weight (with batteries): Approx. 153 grams (5.4 oz)
Power source	2 AA 1.5V alkaline batteries
Pulses, output	Pulse Duration: less than 2 ms Frequency of outputs: 15 to 121 Hz Maximum Voltage: less than 42 V at 500 Ω Maximum Current: less than 500 µA
Waveform	Pulsed current, damped, biphasic, asymmetric sinusoidal with microcurrent output

ADDITIONAL PRODUCT INFORMATION

Modes	Face Stimulate, Face Smooth, Eye Stimulate, Eye Smooth, 60-minute Timeout
User Controls	ON/OFF, Power output intensity increase (\blacktriangleright) and decrease (\blacktriangleleft), Mode selection
Indicators: audible, visual	Mode selection: LED indicates mode selection. Power intensity: LEDs indicate power level. Time: Chirps an beeps. A quiet, subtle tick sound occurs every 30 seconds. A double beep occurs every two minutes. Low battery condition: Indicated by a long series of twenty beeps, in concert with blinking mode LEDs, and followed by device shut off.
Available accessories SOLD SEPARATELY	Y-electrode, brush electrode, pencil electrode, lead wires
Availability, clearances	US FDA cleared for over-the-counter aesthetic use including facial and neck stimulation or body skin stimulation. Cosmetic device available for general purchase. CE approved, CB certified, Health Canada certified. Visit www.avazzia.com for information.
Quality manufacturing	Assembled by Avazzia in U.S. with imported and domestic components. ISO 13485 international certification

WARNING: NOT TO BE USED BY THOSE WITH PACEMAKER OR OTHER IMPLANTABLE ELECTRICAL DEVICE. FOR COMPLETE SAFETY AND USE INSTRUCTIONS, SEE THE OWNER'S MANUAL INCLUDED WITH THE DEVICE.

